WHO TobLabNet Official Method SOP 11

# STANDARD OPERATING PROCEDURE FOR DETERMINATION OF NICOTINE, GLYCEROL AND PROPYLENE GLYCOL IN E-LIQUIDS

No Tobacco Unit (Tobacco Free Initiative)
Tobacco Laboratory Network (TobLabNet)





WHO TobLabNet
Official Method
SOP 11

Standard operating procedure for determination of nicotine, glycerol and propylene glycol in e-liquids

Standard operating procedure for determination of nicotine, glycerol and propylene glycol in e-liquids. WHO TobLabNet Official Method SOP11

ISBN 978-92-4-002274-4 (electronic version) ISBN 978-92-4-002275-1 (print version)

**World Health** 

### © World Health Organization 2021

Some rights reserved. This work is available under the Creative Commons Attribution-NonCommercial-ShareAlike 3.0 IGO licence (CC BY-NC-SA 3.0 IGO; https://creativecommons.org/licenses/by-nc-sa/3.0/igo).

Under the terms of this licence, you may copy, redistribute and adapt the work for non-commercial purposes, provided the work is appropriately cited, as indicated below. In any use of this work, there should be no suggestion that WHO endorses any specific organization, products or services. The use of the WHO logo is not permitted. If you adapt the work, then you must license your work under the same or equivalent Creative Commons licence. If you create a translation of this work, you should add the following disclaimer along with the suggested citation: "This translation was not created by the World Health Organization (WHO). WHO is not responsible for the content or accuracy of this translation. The original English edition shall be the binding and authentic edition".

Any mediation relating to disputes arising under the licence shall be conducted in accordance with the mediation rules of the World Intellectual Property Organization (http://www.wipo.int/amc/en/mediation/rules).

Suggested citation. Standard operating procedure for determination of nicotine, glycerol and propylene glycol in e-liquids. WHO TobLabNet Official Method SOP11. Geneva: World Health Organization; 2021. Licence: CC BY-NC-SA 3.0 IGO

Cataloguing-in-Publication (CIP) data. CIP data are available at http://apps.who.int/iris.

Sales, rights and licensing. To purchase WHO publications, see http://apps.who.int/bookorders. To submit requests for commercial use and queries on rights and licensing, see http://www.who.int/about/licensing.

**Third-party materials.** If you wish to reuse material from this work that is attributed to a third party, such as tables, figures or images, it is your responsibility to determine whether permission is needed for that reuse and to obtain permission from the copyright holder. The risk of claims resulting from infringement of any third-party-owned component in the work rests solely with the user.

**General disclaimers.** The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of WHO concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries. Dotted and dashed lines on maps represent approximate border lines for which there may not yet be full agreement.

The mention of specific companies or of certain manufacturers' products does not imply that they are endorsed or recommended by WHO in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

All reasonable precautions have been taken by WHO to verify the information contained in this publication. However, the published material is being distributed without warranty of any kind, either expressed or implied. The responsibility for the interpretation and use of the material lies with the reader. In no event shall WHO be liable for damages arising from its use.



**No.:** SOP 11

**Date:** 31 March 2021



# World Health Organization Tobacco Laboratory Network

# Standard operating procedure for method

Determination of nicotine, glycerol and propylene glycol in e-liquids

Method: Determination of nicotine, glycerol and propylene glycol

in e-liquids

Analytes: Nicotine (3-[(2S)-1-methylpyrrolidin-2-yl]pyridine)

(CAS # 54-11-5)

Glycerol (propane-1,2,3-triol) (CAS # 56-81-5) Propylene glycol (propane-1,2-diol) (CAS # 57-55-6)

Matrix: e-liquid

Last update: March 2021





**No.:** SOP 11

**Date:** 31 March 2021

### **FOREWORD**

This document was prepared by members of the World Health Organization (WHO) Tobacco Laboratory Network (TobLabNet) in cooperation with member laboratories of the European Joint Action on Tobacco Control (JATC) as an analytical method standard operating procedure (SOP) for measuring nicotine, glycerol and propylene glycol in e-liquids.

### INTRODUCTION

In order to establish comparable measurements for testing e-liquids globally, consensus methods are required for measuring specific contents of e-liquids. The Conference of the Parties (COP) to the WHO Framework Convention on Tobacco Control (WHO FCTC) at its sixth session (Moscow, Russian Federation, 13–18 October 2014) requested the Convention Secretariat to invite WHO to: (a) prepare an expert report on electronic nicotine delivery systems (ENDS) and electronic non-nicotine delivery systems (ENNDS) for the seventh session of the COP (COP7), with an update on the evidence of the health impacts of ENDS/ENNDS, their potential role in quitting tobacco usage and impact on tobacco control efforts; (b) subsequently assess policy options to achieve the objectives outlined in paragraph 2 of decision FCTC/COP6(9); and (c) consider the methods to measure the contents and emissions of these products.<sup>1</sup>

As nicotine content is limited to a certain concentration in some parts of the world (for example, in the European Union, the maximum nicotine concentration in e-liquids is 20 mg/mL), nicotine is considered a priority component to be measured in e-liquids. Since glycerol and propylene glycol are typically ingredients of e-liquids and can be measured simultaneously with nicotine, these components are included in the SOP.

This SOP was prepared to describe the procedure for the determination of nicotine, glycerol and propylene glycol in e-liquid and based on ISO 20714 [2.1].

## 1 SCOPE

This method is suitable for the quantitative determination of nicotine, glycerol and propylene glycol in e-liquids by gas chromatography (GC). The working range of the method is for 1 to 30 mg/mL nicotine, for 200 mg/mL to 1000 mg/mL propylene glycol and for 200 mg/mL to 1000 mg/mL glycerol.

<sup>&</sup>lt;sup>1</sup> Decision FCTC/COP6(9).



### 2 REFERENCES

- 2.1 ISO 20714 (en). E-liquid Determination of nicotine, propylene glycol and glycerol in liquids used in electronic nicotine delivery devices Gas chromatographic method (ISO 20714:2019, IDT).
- 2.2 ISO 13276: Tobacco and tobacco products—Determination of nicotine purity—Gravimetric method using tungstosilicic acid.
- 2.3 ISO 5725-2: Accuracy (trueness and precision) or measurement methods and results—Part 2: Basic method for the determination of repeatability and reproducibility of a standard measurement method.
- 2.4 World Health Organization. Standard operating procedure for validation of analytical methods of tobacco product contents and emissions. Geneva, Tobacco Laboratory Network, 2017 (WHO TobLabNet SOP 02) (https://www.who.int/tobacco/publications/prod\_regulation/standard-operation-validation-02/en/, accessed 10 December 2020).
- 2.5 United Nations Office on Drugs and Crime. Guidelines on representative drug sampling. Vienna, Laboratory and Scientific Section, 2009 (http://www.unodc.org/documents/scientific/Drug\_Sampling.pdf, accessed 10 December 2020).

### 3 TERMS AND DEFINITIONS

- **3.1** *Nicotine content*: total amount of nicotine in e-liquid, expressed as mg per millilitre of e-liquid.
- **3.2** *E-liquid*: liquid or gel which may or may not contain nicotine intended for aerosolization, to be inhaled with an electronic delivery device.
- 3.3 Electronic nicotine delivery device system/electronic non-nicotine delivery system: device used to aerosolize an e-liquid for inhalation.
- 3.4 Laboratory sample: sample intended for testing in a laboratory, consisting of a single type of product delivered to the laboratory at one time or within a specified period.
- 3.5 Test sample: product to be tested, taken at random from the laboratory sample. The number of products taken shall be representative of the laboratory

预览已结束,完整报告链接和二维码如下:

https://www.yunbaogao.cn/report/index/report?reportId=5 23317



