

SAFE BLOOD

AND BLOOD

PRODUCTS

**Manual on
the management,
maintenance and use
of blood cold chain
equipment**



World Health Organization
Geneva

SAFE BLOOD

AND BLOOD

PRODUCTS

Manual on the management, maintenance and use of blood cold chain equipment



World Health Organization
Geneva

WHO Library Cataloguing-in-Publication Data

World Health Organization.

Manual on the management, maintenance and use of blood cold chain equipment.

At head of title : Safe blood and blood products.

1.Blood preservation - instrumentation 2.Blood preservation - methods 3.Plasma 4.Blood platelets
5.Refrigeration - methods 6.Equipment and supplies - standards 7.Guidelines 8.Manuals I.Title

ISBN 92 4 154673 5

(NLM classification: WH 460)

Acknowledgements

The Department of Essential Health Technologies acknowledges the continued support of the Government of Luxembourg towards the WHO Blood Cold Chain Project. The publication was produced by Mr David Mvere, WHO Consultant and Dr Elizabeth Vinelli, National Blood Council, Tegucigalpa, Honduras. Edited by Ms Kay Bond, EHT/WHO. The contribution of the following persons is particularly acknowledged:

Mrs Beryl Armstrong, South African National Blood Service, Pinetown, South Africa
Dr Neelam Dhingra, World Health Organization Headquarters, Geneva, Switzerland
Dr Jean Emmanuel, National Blood Transfusion Service, Malawi
Ms Jan Fordham, World Health Organization Headquarters, Geneva, Switzerland

Reviewers included:

Dr Noryati Abu Amin, World Health Organization Headquarters, Geneva, Switzerland
Dr Rama Bhasin, All India Institute of Medical Sciences, New Delhi, India
Professor Kamel Boukef, National Blood Transfusion Centre, Tunis, Tunisia
Dr Jose Cruz, World Health Organization Regional Office for the Americas/Pan American Health Organization, Washington, D.C., United States of America
Dr Graham Harrison, World Health Organization Regional Office for the Western Pacific, Manila, Philippines
Mrs Shân Lloyd, World Health Organization Headquarters, Geneva, Switzerland
Dr Nishi Madan, University College of Medical Sciences, New Delhi, India
Dr Ana del Pozo, Argentine Association of Hemotherapy and Immunohematology, Buenos Aires, Argentina

© World Health Organization 2005

All rights reserved.

The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of the World Health Organization concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries.

The mention of specific organizations, companies or products does not imply that they are endorsed or recommended by the World Health Organization in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

The World Health Organization does not warrant that the information contained in this publication is complete and correct and shall not be liable for any damages incurred as a result of its use.

Publications of the World Health Organization can be obtained from Marketing and Dissemination, World Health Organization, 20 Avenue Appia, 1211 Geneva 27, Switzerland (Fax: +41 22 791 4857; e-mail: bookorders@who.int). Requests for permission to reproduce or translate WHO publications – whether for sale or for noncommercial distribution – should be addressed to Publications at the above address (fax: +41 22 791 4806; e-mail: permissions@who.int). Information on procuring medical equipment may be obtained from the Department of Contracting and Procurement Services, World Health Organization, 1211 Geneva 27, Switzerland. E-mail: procurement@who.int

Photographs by Sail Vega unless specified

Designed by minimum graphics

Printed in India

Contents

Preface	vii
Useful abbreviations	ix
Glossary	x
1 Introduction	1
1.1 Blood: the raw material	1
1.2 Links in the cold chain	1
1.3 Target audience for this Manual	2
1.4 Using the Manual	3
1.5 Limitations of this Manual	3
2 Storage and transportation of blood and blood components	5
2.1 Safe storage of blood	5
2.1.1 Whole blood	5
2.1.2 Fresh frozen plasma	6
2.1.3 Cryoprecipitate	6
2.1.4 Platelet concentrates	7
2.1.5 Plasma derivatives	8
2.1.6 Cold chain samples and reagents	9
2.2 Packing and transportation of blood and blood components	9
2.2.1 Transportation of whole blood from the collection site to the laboratory	9
2.2.2 Transportation of blood components from one blood bank to another	11
2.2.3 Issuing blood components to clinical areas	12
3 Blood storage equipment: Refrigerators, plasma freezers and platelet agitators	15
3.1 Technical terms for specifications of blood cold chain equipment	15
3.2 Design features common to refrigerators and freezers	17
3.2.1 The cabinet	17
3.3 Ideal design features specific to blood bank refrigerators	19

3.4	Ideal design features specific to plasma (and cryoprecipitate) freezers	21
3.5	Walk-in cold rooms and freezer rooms	21
3.6	Platelet agitators	22
3.7	The cooling mechanism and its maintenance (the refrigeration cycle)	23
3.7.1	Refrigerant gas	24
3.7.2	The compressor	24
3.7.3	The condenser	24
3.7.4	The evaporator	24
3.7.5	The thermostat	25
3.8	Ensuring electrical safety of the equipment	25
3.9	Care of refrigeration equipment	25
4	Other blood cold chain devices	27
4.1	Plasma thawing equipment	27
4.2	Equipment for the transportation of blood	28
4.2.1	Technical terms used	28
4.2.2	Blood transport boxes	28
4.2.3	Ice packs, cooling plates and cooling pouches	30
4.3	Temperature monitoring devices	30
4.3.1	Portable thermometers	31
4.3.2	Temperature recorders/thermographs	31
4.4	Manual recording of temperatures	33
4.5	Alarm systems	34
4.6	Blood warmers	35
4.7	Essential accessories for all refrigeration equipment	35
4.7.1	Voltage regulators (stabilisers)	35
4.7.2	Stand-by generators	35
5	Installing blood refrigerators and plasma freezers	38
5.1	Action on reception of equipment	38
5.2	Siting of refrigerators and plasma freezers	39
5.2.1	Heat and light	39
5.2.2	Air circulation	40
5.2.3	Levelling	40
5.3	Door seals	40
5.4	Cleaning	41
5.5	Energy supply	41
5.6	Starting the equipment	41
5.7	Verifying installations and operational performance	44

6 Organizing the blood cold chain	46
6.1 The structure of a national blood transfusion service	46
6.2 Activities of the blood bank	47
6.3. Critical stages in the movement of blood from collection to transfusion	48
6.3.1 Packing procedures for transportation	49
6.3.2 Receipt and handling on incoming, unprocessed blood and plasma derivatives	50
6.3.3 Receipt and handling of processed blood and blood components	51
6.3.4 Quarantine policies and procedures	51
6.3.5 Labelling of products	52
6.3.6 Method of storage of blood components in available stock	52
6.3.7 Release of whole blood/red cells for use from available stock	53
6.3.8 Procedures for thawing and releasing frozen plasma and cryoprecipitate	53
6.3.9 Procedures for the release of platelet concentrates	53
6.3.10 Discarded blood products and their safe disposal	54
6.4 Monitoring the blood inventory	55
6.4.1 Theoretical count	55
6.4.2 Physical count	56
6.4.3 Daily blood bank report	56
6.4.4 Unused blood components	56
6.5 Model list of essential blood cold chain equipment	59
6.6 Ensuring the blood cold chain during the issuing of blood	60
6.7 Withdrawal of blood from the blood bank, transfusion service or a satellite refrigerator	61
7 Preventive maintenance, care and repair of equipment	63
7.1 Organizing an equipment maintenance programme	63
7.2 Basic care and preventive maintenance of blood cold chain equipment and accessories	67
7.2.1 Blood bank refrigerators and freezers	67
7.2.2 Blood transport boxes	70
7.2.3 Platelet agitators	70
7.2.4 Plasma thawers	70
7.2.5 Stand-by generators	70
7.2.6 Basic preventive maintenance and repair tool kit	71
7.2.7 Calibration of cold chain devices and equipment	72
7.3 Disposal or decommissioning of cold chain equipment	72
8 Monitoring and evaluating the blood cold chain	74
8.1 Definition of terms	74
8.2 Planning for monitoring and evaluation	75
8.3 Quality indicators for evaluation	75

8.4	Records	76
8.5	Reports	76
8.5.1	Blood cold chain performance report	77
8.5.2	Blood cold chain care and preventive maintenance log and repair record	77
8.6	Analysis and interpretation of data	78
8.7	Using the findings from the monitoring and evaluation exercise	78
9	Guidelines for the development of a training programme	80
9.1	Objectives	80
9.2	Key points in the training programme	80
9.3	Preparing a training programme	80
9.4	Developing a training guide	81
	References	84
	Annex 1. The use of CFC in blood cold equipment	85
	Annex 2. WHO Minimum performance specifications for blood cold chain equipment	87
	Annex 3. Basic operational framework for the blood cold chain	92

Preface

The critical contribution that effective management and use of medical equipment brings to health service delivery is only recently gaining recognition. Managing medical equipment has often been misunderstood as the mere procurement of accessible products within a given budget. However, this narrow perspective has proven neither effective, nor cost-effective in the running of health services. The World Health Organization (WHO) promotes the adoption in countries of a comprehensive life cycle approach that falls largely into the following stages: (i) planning and decision-making (e.g. policy, needs assessment and budgeting); (ii) acquisition (including selection, procurement and donation guidelines); (iii) installation (inspection, testing, acceptance, inventories and documentation); (iv) monitoring of performance and use, including preventive maintenance, care and repair; and (v) decommissioning.

The first two stages have been covered extensively in the publication *The Blood Cold Chain: Guide to the Selection and Procurement of Equipment and Accessories*.¹ In that Guide, WHO provided blood bank managers, procurement agencies and manufacturers with a description of, and minimum performance specifications for all the essential equipment needed for the efficient storage and transportation of blood and blood components. WHO plans to update this Guide in line with the improved pre-qualification process for immunization cold chain equipment.²

This new, complementary publication concentrates on the later stages of the life span. Detailed explanations, illustrations and standard operating procedures provide hospital administrators, managers, technicians and all users of blood cold chain equipment with information on how to receive, install, operate, maintain and monitor the equipment. Inspired by

预览已结束，完整报告链接和二维码如下：

https://www.yunbaogao.cn/report/index/report?reportId=5_29922

