

# **Determination of airborne fibre number concentrations**

A recommended method,  
by phase-contrast optical microscopy  
(membrane filter method)



World Health Organization  
Geneva  
1997

WHO Library Cataloguing in Publication Data

Determination of airborne fibre number concentrations : a recommended method, by phase-contrast optical microscopy (membrane filter method).

1. Air pollutants, Occupational — analysis 2. Environmental monitoring — methods  
3. Microscopy, Phase-contrast

ISBN 92 4 154496 1 (NLM Classification: WA 450)

The World Health Organization welcomes requests for permission to reproduce or translate its publications, in part or in full. Applications and enquiries should be addressed to the Office of Publications, World Health Organization, Geneva, Switzerland, which will be glad to provide the latest information on any changes made to the text, plans for new editions, and reprints and translations already available.

© **World Health Organization 1997**

Publications of the World Health Organization enjoy copyright protection in accordance with the provisions of Protocol 2 of the Universal Copyright Convention. All rights reserved.

The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of the Secretariat of the World Health Organization concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries.

The mention of specific companies or of certain manufacturers' products does not imply that they are endorsed or recommended by the World Health Organization in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

Typeset in Hong Kong  
Printed in England

95/10760—Best-set/Clays—7000

# Contents

Preface	v
Outline of the method specification	vii
<b>1. Introduction</b>	<b>1</b>
<b>2. Scope of application</b>	<b>3</b>
<b>3. Specifications of parameters</b>	<b>5</b>
3.1 Sampling	5
3.1.1 Filter	5
3.1.2 Filter holder	7
3.1.3 Storage and transport	9
3.1.4 Sampling pump	9
3.1.5 Optimal filter fibre loadings	10
3.1.6 Flow rate	11
3.1.7 Single sample duration	14
3.1.8 Blanks	15
3.2 Evaluation	16
3.2.1 Sample preparation	16
3.2.2 Optical requirements	22
Microscope	22
Eyepiece graticule	23
Stage micrometer	25
Test slide	26
3.2.3 Counting and sizing fibres	27
Low-power scanning	27
Graticule field selection	28
Laboratory working conditions	28
Counting rules	29
3.2.4 Calculation of fibre concentration	33
<b>4. Accuracy, precision and lower limit of measurement</b>	<b>36</b>
4.1 Accuracy	36
4.2 Precision	37

DETERMINATION OF AIRBORNE FIBRE NUMBER CONCENTRATIONS

4.3 Lower limit of measurement	38
<b>5. Quality assurance</b>	40
References	42
Related reading	44
Annex 1. Static monitoring	45
Annex 2. Characterization of fibres	47
Annex 3. List of participants at the final meeting	52

# Preface

The inhalation of airborne fibres in the workplace can cause a variety of occupational respiratory diseases, which contribute appreciably to morbidity and mortality among workers in both developing and developed countries. Monitoring airborne fibre concentrations is an important tool for occupational health professionals for assessing exposure and establishing the need to control it, evaluating the efficiency of control systems, and characterizing exposure in epidemiological studies. However, a proliferation of methods has hitherto hindered the comparability of results, as well as the possibility of having worldwide proficiency testing to ensure the reliability of obtained data.

A methodological framework that allows for meaningful comparisons between results obtained by different researchers and laboratories is of immense benefit for all aspects of occupational hygiene, but particularly for exposure assessment and environmental monitoring. In addition to proficiency testing, quality assurance schemes, comparisons and exchanges of data and international collaborative studies depend on the use of compatible methodologies. Furthermore, to ensure that preventive control systems in the workplace are adequate and effective, the reliability and comparability of monitoring and exposure data are essential. Occupational health surveillance and the ability to establish correlations between epidemiological and environmental indicators also depend on the ability to compare data from diverse sources.

A project aiming to establish a unified methodology for the evaluation of airborne fibres in the work environment was therefore carried out by the World Health Organization (WHO). A draft text of the present method specification was initially prepared by Dr N. P. Crawford, Institute for Occupational Medicine, Edinburgh, Scotland. During the course of four meetings of an international working group of experts (see Annex 3 for a list of participants at the final meeting),

current evaluation methods were compared and their differences identified and analysed, with a view to understanding the effect of these differences on the results of counting airborne fibres. Consensus was reached by the working group on a recommended method for the determination of airborne fibre number concentrations by phase-contrast optical microscopy. Dr Crawford's role as rapporteur at the final meeting and his work in revising the draft manuscript are gratefully acknowledged.

Collaboration from the International Labour Office (ILO), the European Commission (EC), the International Organization for Standardization (ISO), the European Committee for Standardization (CEN), the International Fibre Safety Group (IFSG), as well as national agencies, particularly the Health and Safety Executive (HSE), United Kingdom, and the National Institute for Occupational Safety and Health (NIOSH), USA, has been fundamental to this project and is also gratefully acknowledged. Special thanks are offered for the financial support provided by the EC and the IFSG.

Particular acknowledgement should also be made of the valuable personal contribution of several members of the working group, particularly Dr N. G. West (HSE), Dr P. Baron (NIOSH), and Mr S. Houston (IFSG), as well as Mrs B. Goelzer, Occupational Health, WHO.

In its first phase, this project has focused on methodology, so that authoritative scientific knowledge can be utilized to ensure accurate and precise measurements of airborne fibre number concentrations. In its second phase, the project will consider worldwide efforts for proficiency testing, quality assurance, and technical cooperation, including training and education.

Dr M. I. Mikheev  
Chief, Occupational Health  
World Health Organization

# Outline of the method specification

## Principle of the method

A sample is collected by drawing a known volume of air through a membrane filter by means of a sampling pump. The filter is rendered transparent (“cleared”) and mounted on a microscope slide. Fibres on a measured area of the filter are counted visually using phase-contrast optical microscopy (PCOM), and the number concentration of fibres in the volume of air is calculated.

## Sampling

Filter:	Membrane of mixed esters of cellulose or cellulose nitrate, 0.8–1.2 $\mu\text{m}$ pore size, 25 mm diameter.
Filter holder:	Fitted with an electrically conducting cowl.
Transport:	In closed holders.
Flow rate:	0.5–16 litres $\cdot\text{min}^{-1}$ . Adjust to give 100–650 fibres $\cdot\text{mm}^{-2}$ .
Blanks:	Sampling media, 4% of filters. Field, $\geq 2\%$ of samples. Laboratory, optional.

## Sample preparation

Acetone–triacetin for fibres with a refractive index  $>1.51$ ; stable for  $\geq 1$  year.

Acetone/etch/water for fibres with a refractive index  $\leq 1.51$ ; unstable.

## Sample evaluation

Technique:	Phase-contrast optical microscopy.
Microscope:	Positive phase contrast, $\times 40$ objective, $\times 400$ –600 magnification.

DETERMINATION OF AIRBORNE FIBRE NUMBER CONCENTRATIONS

	Walton–Beckett graticule, type G-22 ( $100 \pm 2\mu\text{m}$ diameter).
	HSE/NPL Mark II test slide.
	Stage micrometers (1 mm long, $2\text{-}\mu\text{m}$ divisions).
Calibration:	To meet visibility requirements of the test slide.
Analyte:	Fibres (visual count).
Counting rules:	Select counting fields at random, subject to defined criteria.  A countable fibre is $>5\mu\text{m}$ long, $<3\mu\text{m}$ wide and with a length:width ratio $>3:1$ , subject to defined rules when it overlaps the graticule perimeter and when it touches other fibres or particles.
Lower limit of measurement:	10 fibres per 100 graticule fields.
Bias and reproducibility:	See sections 4 and 5.

预览已结束，完整报告链接和二维码如下：

[https://www.yunbaogao.cn/report/index/report?reportId=5\\_30631](https://www.yunbaogao.cn/report/index/report?reportId=5_30631)

