



Rapid Communication:
Molecular assays as initial tests for the diagnosis of
tuberculosis and rifampicin resistance

January 2020

Molecular assays intended as initial tests for the diagnosis of pulmonary and extrapulmonary TB and rifampicin resistance in adults and children: rapid communication. Policy update

ISBN 978-92-4-000033-9

© World Health Organization 2020

Some rights reserved. This work is available under the Creative Commons Attribution-NonCommercial-ShareAlike 3.0 IGO licence (CC BY-NC-SA 3.0 IGO; <https://creativecommons.org/licenses/by-nc-sa/3.0/igo>).

Under the terms of this licence, you may copy, redistribute and adapt the work for non-commercial purposes, provided the work is appropriately cited, as indicated below. In any use of this work, there should be no suggestion that WHO endorses any specific organization, products or services. The use of the WHO logo is not permitted. If you adapt the work, then you must license your work under the same or equivalent Creative Commons licence. If you create a translation of this work, you should add the following disclaimer along with the suggested citation: “This translation was not created by the World Health Organization (WHO). WHO is not responsible for the content or accuracy of this translation. The original English edition shall be the binding and authentic edition”.

Any mediation relating to disputes arising under the licence shall be conducted in accordance with the mediation rules of the World Intellectual Property Organization.

Suggested citation. Molecular assays intended as initial tests for the diagnosis of pulmonary and extrapulmonary TB and rifampicin resistance in adults and children: rapid communication. Policy update. Geneva: World Health Organization; 2020. Licence: [CC BY-NC-SA 3.0 IGO](https://creativecommons.org/licenses/by-nc-sa/3.0/igo).

Cataloguing-in-Publication (CIP) data. CIP data are available at <http://apps.who.int/iris>.

Sales, rights and licensing. To purchase WHO publications, see <http://apps.who.int/bookorders>. To submit requests for commercial use and queries on rights and licensing, see <http://www.who.int/about/licensing>.

Third-party materials. If you wish to reuse material from this work that is attributed to a third party, such as tables, figures or images, it is your responsibility to determine whether permission is needed for that reuse and to obtain permission from the copyright holder. The risk of claims resulting from infringement of any third-party-owned component in the work rests solely with the user.

General disclaimers. The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of WHO concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries. Dotted and dashed lines on maps represent approximate border lines for which there may not yet be full agreement.

The mention of specific companies or of certain manufacturers’ products does not imply that they are endorsed or recommended by WHO in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

All reasonable precautions have been taken by WHO to verify the information contained in this publication. However, the published material is being distributed without warranty of any kind, either expressed or implied. The responsibility for the interpretation and use of the material lies with the reader. In no event shall WHO be liable for damages arising from its use.

Background

Tuberculosis (TB) remains a threat to global public health and is the top infectious cause of death globally. In 2018, an estimated 10 million people developed TB and 1.5 million died from the disease.¹ At least 1 million children become ill with TB every year. In 2018, an estimated 205 000 children died of TB. About 500,000 new cases of multidrug² and rifampicin-resistant tuberculosis (MDR/RR-TB) are estimated to emerge annually but only one in three cases were reported by countries to have been diagnosed and treated in 2018.

Ending the global TB epidemic will only be achievable if there is intensive action by all countries that have endorsed the End TB Strategy and its ambitious targets.³ It requires a paradigm shift from the conventional approaches to disease control used previously. Early and accurate diagnosis of all forms of TB and rapid detection of drug resistance is fundamental to this shift.

Since the approval by WHO of Xpert® MTB/RIF (Cepheid, Sunnyvale, USA, hereafter referred to as 'Xpert MTB/RIF') in 2010, significant additional evidence has been generated on its use as initial test for the diagnosis of TB and rifampicin-resistant TB. In recent months new data on the use of Xpert® MTB/RIF Ultra (Cepheid, Sunnyvale, USA, hereafter referred to as 'Xpert Ultra') and on the latest version Truenat® MTB and MTB Plus system (Molbio Diagnostics, Goa, India, hereafter referred to as 'Truenat') have also become available.

WHO commissioned a systematic review of all the available data in 2019. The results were assessed during a meeting of an independent WHO-convened Guideline Development Group (GDG) on 3-6 December 2019. Detailed recommendations will be published as part of updated WHO Consolidated Guidelines on TB diagnosis in 2020.

This Rapid Communication aims to inform national TB programmes and other stakeholders about the key implications of the latest evidence on the use of specific molecular assays as initial diagnostic tests of pulmonary and extrapulmonary TB and RR-TB, in adults and children.

The updated 2020 WHO Consolidated Guidelines will also incorporate recent WHO recommendations on other rapid tests including line probe assays, urine lipoarabinomannan lateral flow assays and molecular loop-mediated isothermal amplification assays.

Key findings

- High diagnostic accuracy and improved patient outcomes confirmed for Xpert MTB/RIF as initial test to diagnose pulmonary TB (i.e. replacing smear microscopy) and to simultaneously detect rifampicin resistance

¹ Global Tuberculosis Report 2019 https://www.who.int/tb/publications/global_report/en/

² Defined as combined resistance to rifampicin and isoniazid, the two most important anti-TB medicines.

³ Implementing The End TB Strategy: the essentials. Geneva: World Health Organization; 2015 (WHO/HTM/TB/2015.31; http://www.who.int/tb/publications/2015/end_tb_essential.pdf?ua=1)

Intervention: *Xpert MTB/RIF as replacement for sputum smear microscopy (bacteriological culture as reference standard)*

Data assessed: 70 studies (including 5 randomized controlled trials) involving more than 30,000 participants from 37 countries (including 14 high-burden TB countries).

Results: High diagnostic accuracy of Xpert MTB/RIF confirmed in adults with pulmonary TB: Overall sensitivity was 85% in all specimens (smear-positive and smear-negative). Overall specificity was 98% in all specimens (smear-positive and smear-negative). Sensitivity in sputum specimens from HIV co-infected participants was 81% while specificity was 98%.

Results also showed improved patient-important outcomes (cure rates, reduced mortality and reduced pre-treatment loss) when Xpert MTB/RIF replaced microscopy as the initial diagnostic test.

Intervention: *Xpert MTB/RIF for simultaneous detection of rifampicin resistance (phenotypic drug susceptibility testing as reference standard)*

Data assessed: 48 studies involving more than 8,000 participants from 33 countries (including 9 high-burden TB countries).

Results: High overall sensitivity (96%) and specificity (98%) when compared to phenotypic drug susceptibility testing.

- High diagnostic accuracy of Xpert Ultra as initial test to diagnose pulmonary TB (i.e. replacing smear microscopy) and to simultaneously detect rifampicin resistance

Intervention: *Xpert Ultra as replacement for sputum smear microscopy (bacteriological culture as reference standard)*

Data assessed: 6 studies involving more than 2,000 participants from 14 countries (including 8 high-burden countries). More data are needed on relevance of 'trace calls' and on patient outcomes.

Results: High diagnostic accuracy of Xpert Ultra in adults with pulmonary TB: Overall sensitivity (including 'trace' calls as positive) was 90% in all specimens (smear-positive and smear-negative). When 'trace' calls were excluded the overall sensitivity decreased to 86%. Overall specificity (including 'trace' calls) was 96% in all specimens (smear-positive and smear-negative). Excluding 'trace' calls resulted in a slight increase in specificity to 98%. Sensitivity in sputum specimens from HIV co-infected participants was 88% while specificity was retained at 95%.

Intervention: *Xpert Ultra for simultaneous detection of rifampicin resistance (phenotypic drug susceptibility testing as reference standard)*

Data assessed: 5 studies involving close to 1,000 participants from 12 countries (including 7 high-burden TB countries).

Results: High overall sensitivity (94%) and specificity (99%) when compared to phenotypic drug susceptibility testing.

- Use of Xpert MTB/RIF and Xpert Ultra to diagnose TB and detect rifampicin resistance in children from sputum, stool, nasopharyngeal and gastric specimens

Given the difficulties in obtaining sputum specimens from children and the technical limitations of conventional bacteriological methods to aid diagnosis, various non-pulmonary specimens and composite reference standards are often used in evaluating the performance of new diagnostic technologies in paediatric TB.

Intervention: *Xpert MTB/RIF or Xpert Ultra as initial diagnostic test for TB in children (bacteriological culture as reference standard)*

Data assessed for Xpert MTB/RIF: 43 studies involving more than 6,000 participants from 21 countries (including 14 high-burden countries. Studies evaluated sputum specimens, gastric specimens, nasopharyngeal specimen; and stool specimens).

Results: Sensitivity varied by specimen type (46% for nasopharyngeal specimens; 61% for stool; 65% for sputum; and 73% for gastric specimens). Specificity for all specimens varied from 98% to 100%.

Data assessed for Xpert Ultra: Three studies involving close to 700 participants in three high-burden TB countries. These studies all evaluated sputum specimens while one study also evaluated nasopharyngeal specimens.

Results: Sensitivity varied by specimen type (46% for nasopharyngeal specimens; 73% for sputum). Specificity for all specimens varied from 97% to 98%.

Intervention: *Xpert MTB/RIF to detect rifampicin resistance in children (phenotypic drug susceptibility as reference standard)*

Data assessed: Six studies involving more than 200 participants from four high-burden TB countries.

Results: High overall sensitivity (90%) and specificity (98%) when compared to phenotypic drug susceptibility testing.

- Use of Xpert MTB/RIF and Xpert Ultra to diagnose TB and detect rifampicin resistance in adults with extra-pulmonary TB

Given the difficulties in obtaining extra-pulmonary specimens and the technical limitations of conventional bacteriological methods to aid diagnosis, various non-pulmonary specimens and composite reference standards are often used in evaluating the performance of new diagnostic technologies in extra-pulmonary TB.

Intervention: *Xpert MTB/RIF as initial diagnostic test in adults with extra-pulmonary TB (bacteriological culture as reference standard)*

Data assessed: 59 studies from 26 countries (including seven high-burden countries). MTB/RIF in adults with extra-pulmonary TB. Studies evaluated cerebrospinal fluid specimens, lymph node aspirate, lymph node biopsy, pleural fluid, urine, synovial fluid, peritoneal fluid, pericardial fluid, and blood.

Results: Sensitivity of Xpert MTB/RIF varied by specimen type (from 50% for pleural fluid up to 97% for synovial fluid). The specificity of Xpert MTB/RIF also varied by specimen type (from 79% for lymph node biopsy up to 99% for pleural fluid).

Intervention: *Xpert Ultra as initial diagnostic test in adults with extra-pulmonary TB (bacteriological culture as reference standard)*

Data assessed: Six studies involving from three countries (including two high-burden TB countries). Studies evaluated cerebrospinal fluid, pleural fluid, urine, synovial fluid and urine.

Results: Sensitivity of Xpert Ultra varied by specimen type (from 71% for pleural fluid up to 100% for lymph node biopsy). The specificity of Xpert Ultra also varied by specimen type (from 71% for pleural fluid up to 100% for urine).

Intervention: *Xpert MTB/RIF or Xpert Ultra for simultaneous detection of rifampicin resistance (phenotypic drug susceptibility testing as reference standard)*

Data showed overall high performance of both assay for detection of rifampicin resistance: sensitivity 96-97% and specificity 99%.

- High diagnostic accuracy of Truenat as initial test to diagnose TB (i.e. replacing sputum smear microscopy) and to sequentially detect rifampicin resistance

Intervention: *Truenat MTB and MTB Plus as replacement for sputum smear microscopy (bacteriological culture as reference standard)*

Data assessed: Interim data from a FIND-coordinated multi-central, prospective field evaluation study in

预览已结束，完整报告链接和二维码如下：

https://www.yunbaogao.cn/report/index/report?reportId=5_24852

